

## Major strengths

- Easy to use
- Own two internal shelves
- Light warning during heating
- Timer goes from 1 to 120 minutes
- Sized for a mold
- Internal chamber in aluminum for a better heat diffusion
- Adjustable flow air going inside the chamber for venting

## A device adapted for Soft-lithography applications

- Size adapted to size of a Petri dish to heat two molds at one time but with the minimum place and electric consumption.
- It is pre-heated rapidly.
- It is a really movable device and can be install everywhere.

## Safety & protection

- The level protection is IP40.
- The quality of the product is certified at an international level.

## Warranty

Warranty 1 year.



BlackHole·Lab

They trust us



# ▶ SoftLithoBox®

## Oven

▶ **FEATURES**

## Oven

Natural convection drying ovens with notable thermal stability. Frontal panel with bimetallic fluid expansion probe thermo regulator. This oven sized for 4 inches wafer will allow you to cure your PDMS efficiently.



▶ **SPECIFICATIONS**

Weight	10 kg (22 lbs)
Power	360 W
Power supply	110V or 230V; 50/60 Hz
Dimensions (HxWxD)	380x240x300 (in mm)
Internal dimension (HxWxD)	190x150x190 (in mm)
Temperature range	Ambient +1°C to 150°C
Temperature resolution	3°C
Volume	5.4 liter

Phone: 0033(0).782.831.207  
E-mail: [contact@blackholelab.com](mailto:contact@blackholelab.com)  
[www.blackholelab.com](http://www.blackholelab.com)

BlackHole Lab  
111 avenue Victor Hugo  
75784 Paris cedex 16