

Major strengths

- Ultra-quiet
- Long-life
- Easy to operate and minimum maintenance
- Handy and easily portable
- Anti-corrosion treated valve plates
- Best-in-class, long-life synthetic oil

A device adapted for Soft-lithography applications

- With a minimum air flow of 1 l/min and a pressure of 8 bar it will totally fit to supply your spin coater or your air gun.
- Its little size and its little noise enables to install it everywhere without bothering users.
- The add of a air and oil filter protect the wafer of oil residue and particle pollution.

Safety & protection

- A safety valve prevents the pressure to go over 16 bar.
- The vessel can resist at a pressure of 24 bar.

Warranty

Warranty 2 years.



BlackHole·Lab

They trust us



▶ SoftLithoBox®

Air compressor

▶ FEATURES

Air compressor

This lubricated air compressor is very easy to operate. It is little and noiseless and so really discreet. Made as much for industrial application as laboratory, it is a really versatile tool.



▶ SPECIFICATIONS

Weight	18 kg (39.6 lbs)
Power	350 W
Power supply	120V or 230V 50/60 Hz
Dimensions (HxWxD)	380x333x342 (in mm)
Volume tank	4 L
Maximum pressure	8 bar (120 PSI)
Airflow	From 1 l/min to 20 l/min
Noise level	< 35 dBA

Phone: 0033(0).782.831.207

E-mail:
contact@blackholelab.com
www.blackholelab.com

BlackHole Lab
111 avenue Victor Hugo
75784 Paris Cedex 16