

Majors strengths

- Compact tool easily movable.
- Short vibration orbit..
- Adjustable agitation frequency.
- Easy of use.

A device adapted for Soft-lithography applications

- Enables a control and repeatable development of your photoresist.
- Its short vibration orbit enables a homogenous and perfect development and to remove bubbles that can be trap close to the design.

Safety & protection

A large selection of attachments enables to use various vessels and sizes without any risks.

Warranty

Warranty 3 years.



BlackHole·Lab

They trust us



▶ SoftLithoBox®

Developer agitator

▶ FEATURES

Developer agitator system

This tool is a vibrating agitator, electronically regulate it allows smooth to ridged mixing of viscous media and solids.

Already used by different laboratory for microfluidic application, it enables an easy photoresist development.



▶ SPECIFICATIONS

Weight	5.5kg (12.1lbs)
Power	31 W
Power supply	110V or 230V ; 50 Hz
Dimension (HxWxD)	125x245x310 (in mm)
Weight capacity	2 kg
Vibration orbit	3 mm
Speed	From 150 to 1350 rpm
Vibration time	From 1 to 120 minutes
Size plate	220x220 (in mm)

Phone: 0033(0).782.831.207
E-mail: contact@blackholelab.com
www.blackholelab.com

BlackHole Lab
111 avenue Victor Hugo
75784 Paris cedex 16